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
The bottom of the page features several handwritten signatures in black ink. To the right of the signatures is a blue circular stamp with the text 'सचिव' (Secretary) in the center. Further to the right is a red circular stamp with the text 'नेपाल सरकार' (Government of Nepal) at the top, 'कृषि तथा प्रशुध्नी विकास मन्त्रालय' (Ministry of Agriculture and Animal Development) in the middle, and 'सिन्धुदरबार, काठमाडौं' (Sindhudharbar, Kathmandu) at the bottom.

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
The bottom of the page features several handwritten signatures in black ink. To the right of the signatures is a blue circular stamp with the text 'सचिव' (Secretary) in the center. Further to the right is a circular official seal of the Government of Nepal, with the text 'नेपाल सरकार' (Government of Nepal) and 'कृषि तथा पशुपन्छी विकास मन्त्रालय' (Ministry of Agriculture and Animal Development) around the perimeter. Below the seal, the text 'सिंहवजार, काठमाडौं' (Sinhvajar, Kathmandu) is visible.

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


Handwritten signatures and a blue stamp reading "सचिव" (Secretary).

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1. Particle Size Distribution (Mechanical Analysis)

Apparatus

1. Soil hydrometer ASTM graduated (-5 to 60)
2. Hydrometer jar
3. Mechanical stirrer
4. Dispersion cup
5. Beaker 250 ml
6. Pipette, 10 ml

Reagents:

1. Sodium Hexametaphosphate: Dissolve 101 gm. of Sodium Hexametaphosphate in 1 L water.

Procedure:

Weigh **50 gm.** of soil in a 250 ml beaker and add sufficient water to cover the soil. Then add **10 ml** of sodium hexametaphosphate solution, stir well with a glass rod and leave overnight. Transfer it in a dispersion cup and add sufficient water to fill two- third of the cup. Stir for 10 minutes in the mechanical stirrer, transfer in the hydrometer jar and make up the volume to the mark with the hydrometer in it. Remove the hydrometer and shake the jar upside down several times closing the mouth either by hand or cork. When the soil is well dispersed keep it in the table and note the time immediately. Immerse the hydrometer in the jar and read it at 40 sec and **2 hours**.

(Once in the table it should not be disturbed throughout the experiment). Note the temperature of the suspension at the time of taking hydrometer readings. Correct the hydrometer reading by subtracting 0.3 for every °C below 20°C or by adding if the temperature is above 20°C.

Calculation

(Silt + Clay)% = Hydrometer reading at 40 sec + $0.3 \times (t-20)^\circ \text{C}$

Clay % = Reading at 2 hr. + $0.3 \times (t-20)^\circ \text{C}$

%Sand = $100 - \%(\text{Clay} + \text{Silt})$

%Silt = $\%(\text{Clay} + \text{Silt}) - \% \text{Clay}$

Determine the texture from the triangular chart.

Handwritten signatures and stamps are present at the bottom of the page. A blue stamp reads "सचिव" (Secretary). A circular official stamp of the Government of Nepal is also visible, with text in Nepali: "नेपाल सरकार" (Government of Nepal), "श्री तपो पशुपती विकास योजना" (Shri Tapo Pasupati Vikas Yojana), "सिंहदरबार, काठमाडौं" (Sindharbar, Kathmandu).


Figure: Textural Triangle

1. pH DETERMINATION

Standard Buffer:

1. p^H 4.0: Dry potassium biphthalate ($KHC_8H_4O_4$) for 2 hours at $105^\circ C$ and dissolve 10.21 gm. of it in hot distilled water and dilute the solution to 1 L with distilled water. As a preservative, add 1.0 ml of chloroform or toluene. This solution has a p^H value of 4.0 at temperature from $15-30^\circ C$.
2. p^H 6.86: Dry the two salts KH_2PO_4 and Na_2HPO_4 for 2 hr. at $105^\circ C$ and dissolve 3.44 gm. of KH_2PO_4 and 3.55 gm. of Na_2HPO_4 in distilled water and dilute to 1 L. As a preservative add 1.0 ml of chloroform or toluene. This solution has p^H 6.90 at $15^\circ C$, 6.88 at $20^\circ C$ and 6.85 at $30^\circ C$.

Apparatus:

1. p^H meter
2. Beaker 50 ml
3. Mechanical stirrer
4. glass rod

Procedure:

Weigh 10 g of air dried soil samples (< 2 mm) in 50 ml beaker and add 25 ml of distilled water. Shake for about 1 minute in a mechanical stirrer and leave for about an hour. The p^H meter is calibrated by using standard buffer solution of p^H 4.0 and 6.86. The soil water suspension is stirred well with a glass rod just before immersing the electrode. The p^H of

soil water suspension is then measured with the calibrated p^H meter. For the peat and muck soil, volume of water may be increased

Table: Rating of soil according to p^H

Soil Reaction (p^H)	p^H range
extremely acidic	<5
strongly acidic	5.0-5.5
moderately acidic	5.5-6.0
slightly acidic	6.0-6.5
nearly neutral	6.5-7.5
slightly alkaline	7.5-8.0
moderately alkaline	8.0-8.5
strongly alkaline	8.5-9.0
extremely alkaline	>9.0

2. AVAILABLE POTASSIUM

a) flame photometer method:

Apparatus

- 1) flame photometer
- 2) shaking apparatus
- 3) 100 ml. Erlenmeyer flask

Reagents

- 1) 1 N Ammonium acetate pH 7.0: Dissolve 77 gm. of ammonium acetate in 1L of distilled water.

Or

To 58 ml of glacial acetic acid add 500 ml. of distilled water followed by 65 ml of liquor ammonium and dilute to 1L. Adjust the pH to 7.0 ± 0.02 with dilute NH_4OH or acetic acid.

- 2) K standards (stock solution): Dissolve 0.1905 gm. dried KCl in 1L volumetric flask and make up the volume- 100 ppm K

Take 0, 5, 10, 15, 20 and 25 ml of 100 ppm K solution in 100 ml volumetric flask and dilute with 1N ammonium acetate pH 7.0 solution to the mark- 0, 5, 10, 15, 20 and 25 ppm K.

Procedure

Weight 2 gm. of air dried soil in a 125 ml conical flask (or 100 ml polyethylene bottle with leak proof cap), add 20 ml normal neutral ammonium acetate, shake for 5 minutes in a mechanical shaker and filter through Whatman No. 42, 12.5 cm filter paper. Prepare a standard curve of K by aspirating 0, 5, 10, 15, 20 and 25 ppm K after adjusting full scale deflection of flame photometer with 25 ppm K. Note the reading and draw graph. Aspire the soil solution, note its reading and determine K in the soil solution from the graph.

Calculation

$$K_2O \text{ (kg/ha)} = R \times 20/2 \times 1.2 \times 2.24$$

Where,

R= potassium of soil extract in ppm, from the standard curve

1.2= conversion factor for K to K_2O

2.24 = conversion factor for ppm to kg/ha

20/2= dilution factor

Interpretation

Low- <110 kg/ha K_2O

Medium- 110-280 kg/ha K_2O

High- >280 kg/ha K_2O

3. AVAILABLE PHOSPHORUS

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Modified Olsen's bicarbonate method

Principle

This method uses 0.5 N NaHCO_3 solution of pH 8.5 as an extractant which controls the activity of Ca, Al and Fe by precipitating calcium as carbonate and aluminum and iron as hydroxides. The organic matter dissolved by the extractant must be removed by the use of activated charcoal.

Apparatus

1. 100 ml polythene bottles
2. shaker
3. funnel
4. what man no 42 filter paper
5. volumetric flask 50 ml
6. pipettes 5 and 10 ml
7. beaker 50 ml

Reagents

1. Extracting solution (0.5 N NaHCO_3 , pH 8.5): Dissolve 210 gm of CP NaHCO_3 in 5L distilled water. Adjust the p^{H} to 8.5 with 0.5 N NaOH or H_2SO_4 . As the p^{H} of the solution tends to increase on exposed to atmosphere, few drops of liquid paraffin should be added and the p^{H} should be checked monthly.
2. 5N H_2SO_4 : Dilute 35 ml concentrated A.R. H_2SO_4 to 250 ml.
3. Ammonium Molybdate:
 - a) Dissolve 12 gm. A.R. ammonium molybdate in 250 ml distilled water. In 100 ml distilled water dissolve 0.2908 gm. of antimony potassium tartarate. Add both the solutions in 1000 ml of 5N H_2SO_4 (141 ml of concentrated H_2SO_4 per liter water), mix thoroughly and make to 2L. Store in a pyrex glass bottle and keep it in a dark and cool temperature.
 - b) Dissolve 1.056 gm. of ascorbic acid in 200 ml of ammonium molybdate solution (Reagent 3a). This reagent should be prepared as required since it cannot be kept for more than 24 hours.
4. Activated Charcoal (Darco G- 60)
5. Standard P solution
Primary standard 50 ppm P: 0.2195 gm. of A.R. KH_2PO_4 dried at 40°C is dissolved in about 400ml of distilled water in one liter volumetric flask. Add 25 ml of 7N H_2SO_4 to it and make the volume to 1 L. Thus preserved with H_2SO_4 , the solution keeps indefinitely but it should be stored in a weathered soft glass bottle (rather than one of pyrex) to minimize contamination with arsenate.
Secondary standard 5 ppm P: Dilute 5 ml of 50 ppm P stock solution to 50ml in a volumetric flask for the 5 ppm P. It must be made up if fungus growth in the solution is noticed.
6. p-nitrophenol indicator 0.25: Dissolve 0.25 gm. indicator in 100 ml of distilled water.

Procedure

Weigh 2 gm. soil sample (air dried < 2 mm) in a 100 ml polyethylene bottle. Add one teaspoon of activated charcoal and 40 ml of 0.5 N NaHCO_3 , extracting solution. Shake for 30 minute in a shaker and filter through Whatman no 42 filter paper. Pipette 10 ml aliquot of the filtrate in a 50 ml volumetric flask and acidify with 5N H_2SO_4 to p^{H} 5.0 using p-nitrophenol indicators till the yellow color just disappeared. Shake gently after each addition of acid. Further add acid, drop wise this time until the color changes from yellow to colorless. Add distilled water washing down the sides of volumetric flask to 40 ml followed by 8 ml of reagent (ammonium molybdate). Make up the volume to the mark and shake well. Maximum intensity of the blue color is obtained in 10 minutes and remains stable up to 24 hrs. Include a blank in every batch by shaking the extracting

Handwritten signatures and initials.

solution without soil. It should include all the reagents added in every step. Measure color intensity in a colorimeter after 10 minutes using red filter (660 nm). Prepare the standard curve by taking 0, 0.5, 1, 2, 4, 6, 8, and 10 ml of 5 ppm standard solution 50 ml volumetric flask, and NaHCO_3 extracting solution and proceed exactly like in the solution.

Calculation

$$\text{P}_2\text{O}_5 \text{ (kg/ha)} = \text{ppm P in solution} \times 2.29 \times 100 \times 2.24$$

Where,

100 = dilution factor

2.24 = conversion factor for ppm in soil to kg/ha in soil

2.29 = conversion factor for P to P_2O_5

Interpretation

Low	< 31 kg/ha P_2O_5
Medium	31-55 kg/ha P_2O_5
High	> 55 kg/ha P_2O_5

ORGANIC MATTER DETERMINATION (WALKLEY- BLACK METHOD)

Apparatus:

50 ml. Burette

10 ml. bulb pipette

Measuring cylinder 25 ml. or Acid Dispenser or pipette 20 ml. capacity

Conical flask 500 ml. capacity

Reagents for Organic Matter

- 1) sodium fluoride
- 2) 1 N Potassium Dichromate ($\text{K}_2\text{Cr}_2\text{O}_7$): Dry $\text{K}_2\text{Cr}_2\text{O}_7$ in the temperature of 105°C . Weight 49.04 gram of it and dissolve it in distilled water to make the volume of 1 liter.
- 3) 0.5 N Ferrous ammonium sulphate (FAS): Dissolve 196 gm. of FAS in 800 ml of distilled water. Add 20 ml of conc. sulphuric acid and make volume to 1 liter.
- 4) Ferroin indicator: Dissolve 1.5 gm. of 1-10 Phenanthroline monohydrate ($\text{C}_{12}\text{H}_8\text{N}_2\text{H}_2\text{O}$) and 0.7 gm ferrous sulphate heptahydrate ($\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$) in 100 ml of distilled water.

Procedure

Weigh 0.5gm. Soil sample, passing through 0.2 mm sieve and add exactly 10 ml. of 1 N $\text{K}_2\text{Cr}_2\text{O}_7$ solution to it, in a 500 ml. conical flask. Add 20 ml. of concentrated sulphuric acid and mix by gentle rotation for 1 minute, to insure complete contact of the reagent with the soil, but with care to avoid throwing up soil to the sides of the flask. The mixture is allowed to stand for 30 minutes. A standardization blank (without soil) is run in the same way. After half an hour, add about 200 ml of distilled water, 10 drops of ferroin indicator and about 0.2 gm. of sodium fluoride. Back titrate the solution with ferrous ammonium sulphate solution. The color is dull green with chromic ion in the beginning, and then shifts to green as the titration proceeds. At the end point the color sharply shifts to brick red. If over 8 ml of 10 ml. chromic acid has been consumed during the titration, the determination is repeated with a small quantity of soil or by adding double or treble amount of $\text{K}_2\text{Cr}_2\text{O}_7$ and H_2SO_4 .

$$\text{OM}\% = \frac{(B-S) \times N \times 3 \times 100 \times 100 \times 100}{\text{wt of soil} \times 1000 \times 77 \times 58} = \frac{(B-S)N}{\text{wt of soil}} \times 0.67$$

Where,

B = Volume of FAS used up for blank titration

S = Volume of FAS used up for sample titration

N = Normality of FAS from blank titration

Wt. = weight of soil sample

Note:

- i) Organic matter contains 58 percent carbon.

[Handwritten signatures and marks]

- ii) Recovery factor for this method is 77 percent.
- iii) Equivalent weight of carbon is 3.

1) TOTAL NITROGEN

Apparatus:

1. Kjeldhal digestion flask – 50 ml or block digester tubes
2. N- Digestion apparatus
3. Distilling apparatus
4. Volumetric flask – 100 ml.
5. Conical flask – 125 ml.
6. Acid dispenser or measuring cylinder
7. Pipettes, - 10 and 20 ml.
8. Burette – 25 ml.
9. Kjeldhal flask holder
10. Asbestos Glove

Reagents

1. Digestion Mixture (Catalyst): Grind and mix 10 gm. of copper sulphate with 200 gm. of sodium sulphate.
2. Concentrated sulphuric acid – L.R.
3. Sodium hydroxide: Dissolve 400 gms of sodium hydroxide (flakes or pellets L.R.) in one liter of distilled water and cool.
4. Mixed indicator: Dissolve 0.5 gm. bromo-cresol green and 0.1 gm. methyl red in 100 ml. of 95 percent ethanol.
5. Boric acid 4%: Dissolve 40 gm. boric acid crystal (H_3BO_3) in one liter distilled water.
6. 0.01N HCl: Dilute 17 ml. conc. HCl with distilled water to 2L-(A) Standardize 20 ml of solution (A) with Na_2CO_3 and calculate its normality. Dilute the solution (A) according to its strength to give 0.01 N HCl. [or- Dilute 100 ml. of solution (A) to 1L and calculate its strength from the above normality of (A)]

Procedure:

Weigh 1 gm. soil in 50 ml Kjeldhal digestion flask and add 2 gm. catalyst digestion mixture followed by 10 ml conc. H_2SO_4 and few pieces of broken porcelain. For a fine textured soil add 10 ml. of distilled water and leave it for 30 minutes before adding digestion mixture and sulphuric acid. Mix the soil with sulphuric acid by swirling the flask and heat in the low heat until frothing stops. Gradually increase the heat until the acid boils. Swirl the flask at intervals and digest until the color of the mixture changes to green- blue or grey color and continue it for 1- 1.5 hours more. Care should be taken during digestion not to allow the flame to touch the flask above the part occupied by the liquid.

With Block Digester Decator, add 10 ml concentrated sulphuric acid to 1 gm. soil and 2 gm. digestion mixture in 250 ml digestion tube and heat in preheated Block Digester at $410^{\circ}C$ for 45 minutes. If exhaust system is used, adjust its rate to maximum in the beginning of the digestion and reduce it after about 10 minutes (when the evolution of dense fumes of sulphuric acid decrease) so that the acid fume will be condensed at about two- third of the digestion tube.

Cool the flask and add about 20 ml. of distilled water before the solution starts crystallizing. Transfer the solution in a 100 ml. volumetric flask, leaving the sand in the digestion flask and make up the volume. Take 20 ml. aliquot in the distilling flask and add 20 ml. of 40% NaOH and distil it, collecting the liberated NH_3 in 10 ml 4% boric acid solution containing 2 drops of mixed indicator in 125 ml. conical flask. Titrate it with 0.01 N HCl. Run a blank for each batch of 12 samples.

(Handwritten signatures and initials)

Calculation

$$\%N = \frac{7 \times n \times (T - B)}{S}$$

Where, n= Normality of acid

T= vol. of acid used in titration

B= vol. of acid used in blank

S= Sample wt.

Interpretation

a) Terai

Low.....<0.075%N

Medium.....<0.075- 0.150%N

High.....>0.150%N

b) Hills

Low.....<0.1%N

Medium.....<0.1- 0.3%N

High.....>0.3%N

MICRONUTRIENTS

Apparatus:


1. Atomic Absorption Spectrophotometer
2. Acetylene gas cylinder
3. Paraffin wax paper
4. Shaker
5. Pipette 20 ml
6. Conical flask 125 ml

Reagents:

1. Diethylenetriaminepenta acidic acid (DTPA) extracting solution:
Dissolve 149.2 gm of reagent grade triethanolamine (TEA), 19.67 gm. of DTPA and 14.7 gm. of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ in deionized distilled water (Glass double distilled). Since DTPA is not very soluble in water, place it in a small amount of water and then dissolve in the TEA solution.) When the DTPA has dissolved, dilute to approximately 9 liters. Adjust the pH to 7.3 ± 0.05 with 1:1 HCl (approximately 8.3 ml of 1:1 HCl are required). Make to 10L final volume. The pH should be checked periodically, because a pH of 7.3 is critically for the extraction. Thus prepared 0.005 M DTPA, 0.01 M CaCl_2 and 0.1 M TEA solution is stable for several months.
2. Glass double distilled water or deionized distilled water (DDW).
3. Standard Zinc solution:
 - i) Stock Solution: Dissolve exactly 0.4399 gm. AR $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ in glass double distilled water, add 25 ml 1 N HCl and dilute to 1 L volumetric flask – 100 ppm Zn.
 - ii) Working Standard: They are prepared by diluting 50 ml of 100 ppm Zn to 250 ml (20 ppm Zn) and then further dilution of 0, 2.2, 5, 10, 15 and 25 ml of this 20 ppm Zn to 250 ml with DTPA extracting solution to give 0, 0.2, 0.4, 0.8, 1.2 and 2 ppm of Zn.
4. Standard Fe solution
 - i) Stock solution: Dissolve 0.7023 gm. AR $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$ in deionized distilled water (DDW) and add 20 ml of 1 N H_2SO_4 . Oxidize it by adding 25 ml 1% KMnO_4 slowly and then drop wise till the pink color just stays. Dilute to exactly 1L -100 ppm Fe.
 - ii) Working Standards: They are prepared by diluting 0, 5, 10, 15, 20 and 25 ml of 100 ppm Fe to 100 ml with DTPA extracting solution to give 0, 5, 10, 15, 20 and 25 ppm Fe.

Handwritten signatures and initials.

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5. Standard Copper Solution

- i) Stock Solution: Dissolve 0.3929 gm. AR $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in 500 ml DD water, add 20 ml 1N H_2SO_4 and dilute to 1 L in volumetric flask – 100 ppm Cu.
- ii) Working Standards : They are prepared by diluting 50 ml of 100 ppm Cu to 250 ml (20 ppm Cu) and then further dilution of 0, 2.5, 5, 10, 15 and 25 ml of this diluted (20 ppm Cu) to 250 ml with DTPA extracting solution to give 0, 0.2, 0.4, 0.8, 1.2 and 2.0 ppm of Cu.

6. Standard Mn Solution

- i) Stock Solution: Dissolve 0.2877 gm. potassium permanganate in 500ml DD water and add 25 ml 1N H_2SO_4 . Boil for few minutes and carefully add 1.2 gm sodium sulphate crystal ($\text{Na}_2\text{SO}_4 \cdot 7\text{H}_2\text{O}$). Boil again to remove SO_2 , cool and make up the volume to 1L- 100 ppm Mn.
- ii) Working Standards: They are prepared by diluting 0, 5, 10, 15, 20 and 25 ml of 100 ppm Mn solution to 100 ml with DTPA extracting solution to give 0, 5, 10, 15, 20 and 25 ppm Mn.

Extraction Procedure:

Weigh 10 gm. of air dried soil in a 125 ml Erlenmeyer flask and add 20 ml of the DTPA extracting solution. Cover the flask with paraffin or polyethylene stopper and secure upright on a horizontal shaker with a stroke of 8.0 cm and a speed of 120 cycles/ min. After exactly 2 hours of shaking, filter the suspension by gravity through Whatman No.42 filter paper. The shaking time is very important because extracting is not complete in 2 hours. (The labile and non-labile quantities of the trace elements will continue to dissolve even after 2 hours). The filtrate may be used to analyze for the trace elements Fe, Cu, Mn, and Zn with AAS after standardization for particular element.

a) Available Zinc

Place Zinc Hollow Cathode Lamp in the working lamp turret and let it warm for specified period as per the instruction of the manufacturer. Light the burner using air and acetylene gas, adjust the slit, wavelength and burner height as per the instruction of AAS manufacturer. Let it warm for 5 minutes and standardize AAS using 0 and 2 ppm Zn standard. Run series of Zinc standards, note the absorbance reading and draw graph.

Aspirate the soil extract, note the absorbance reading and determine zinc in solution from the graph.

Calculation:

$$\text{ppm Zn in soil} = 2 \times (R-B)$$

Where,

R is the ppm Zn in soil extract from the graph.

B is the ppm Zn in Blank extract from the graph.

Interpretation

Factors affecting the interpretation of the soil test data are pH, CaCO_3 percent, phosphorus, organic matter, clay percent, CEC.

Typical range is critical level for DTPA extractable Zn is reported to be 0.2 to 2.0 Zn mg/kg (ppm) and average critical value is 0.8 Zn mg/kg.

b) Available Copper

Place copper Hollow Cathode Lamp in the working turret, proceed as in the Zinc determination, note the absorbance reading of copper standards and draw graph.

Aspirate the soil extract, note the absorbance reading and determine Cu in solution from the graph.

Calculation

$$\text{ppm Cu in soil} = 2 \times (R-B)$$

Where,

[Handwritten signatures and initials]

R is the ppm Cu in soil extract from the graph.
B is the ppm Cu in Blank extract from the graph.

Interpretation

Factors affecting the interpretation of the soil test data are pH, CaCO_3 percent and crop.

Typical range is critical level for DTPA extractable Cu is reported to be 0.12 to 0.25 Cu mg/kg (ppm).

c) Available Iron

Place copper Hollow Cathode Lamp in the working turret, proceed as in the Zinc determination, note the absorbance reading of Iron standards and draw graph.

Aspirate the soil extract, note the absorbance reading and determine Fe in solution from the graph.

Calculation

ppm Cu in soil = $2 \times (R-B)$

Where,

R is the ppm Fe in soil extract from the graph.

B is the ppm Fe in Blank extract from the graph.

Interpretation

Factors affecting the interpretation of the soil test data are pH, CaCO_3 percent, organic matter, aeration, soil moisture and CEC.

Typical range is critical level for DTPA extractable Fe in USA and Tropics is reported to be 2.5 to 5.0 Fe mg/kg(ppm).

d) Available Manganese

Place copper Hollow Cathode Lamp in the working turret, proceed as in the Zinc determination, note the absorbance reading of copper standards and draw graph.

Aspirate the soil extract, note the absorbance reading and determine Mn in solution from the graph.

Calculation

ppm Cu in soil = $2 \times (R-B)$

Where,

R is the ppm Mn in soil extract from the graph.

B is the ppm Mn in Blank extract from the graph.

Interpretation

Factors affecting the interpretation of the soil test data are pH, organic matter, CaCO_3 percent and texture.

Typical range is critical level for DTPA extractable Mn is reported to be 1.0 to 5.0 Mn mg/kg (ppm) and average critical value is 1.4 mg/kg.

e) Available Boron (hot water extraction method)

Chemical and reagents

1) Buffer

35 g ammonium acetate

+

2.5g potassium acetate

+


1g nitrilotriacetic acid


+

Dissolve in 100ml of DH_2O


+

Slowly add 13.25 ml of glacial acetic acid and make volume to 250ml.

२२

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- 2) **Extraction solution (calcium chloride 0.02N)**
Dissolve 1.11gm of anhydrous CaCl_2 in 500 ml of distilled water and dilute to exact 1000 ml.

- 3) **Azomethine (H)**

0.3 g azomethine H

+

2g ascorbic acid

+

Dissolve in 60 ml DDH_2O and make volume to 100ml.

+

It should be prepared at least 24 hours before use and store in polypropylene bottle wrapped in a foil and place in refrigerator.

- 4) **Standard boron solution - 100ppm**

0.5716 g AR boric acid in 1 liter DDH_2O gives 100 ppm B.

Take 0, 1, 2, 3, 4 and 5 ml of 100 ppm B in 100 ml volumetric flask and make volume by DDH_2O . It gives 0, 1, 2, 3, 4 and 5 ppm B.

Boron analysis procedure

10g soil

+

20 ml extraction solution

+

5 minute heating in boiling hot water

+

Cool in a running tap and filter

+

Cut the sealing of the polypropylene plastic paper

+

Pour the solution to Filter through Whatman 42 filter paper

+

Take 1ml aliquot in PVC tube

+

Add 2ml buffer solution

+

2ml azomethine H reagent

+

Mix well and take reading in spectrophotometer at 420nm after 30 minutes

+

Prepare standard with 0, 1, 2, 3, 4, and 5 ppm B and develop color as above and take readings.

Put in the packet of sealed polypropylene plastic paper

Calculation

ppm B in soil = ppm of B in solution X dilution factor
= ppm of B in solution X (2 X 5)

f) Available Molybdenum

Method: Acid Ammonium Oxalate

Apparatus

- Erlenmeyer flask or narrow mouth polyethylene bottle

[Handwritten signatures]

[Handwritten signature]

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- Shaker
- Separating funnel with rack
- Spectrophotometer
- Water bath
- Muffle furnace

Reagents

- Glass double distilled water or deionized distilled water
- Extracting solution Acid Ammonium Oxalate, 0.2N $\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$ with 2.5 percent $\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$, pH 3.3: Dissolve 24.9g ammonium oxalate and 12.6 g oxalic acid in 1 L distilled water. The pH of the solution should be near 3.3.
- Hydrochloric acid 6.5N, containing Fe^{3+} : Dissolve 0.5g of $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ in 1 Litre of 6.5N HCl
- Citric acid (Crystal)
- Organic extractant : Mix equal volume of carbaon tetra-chloride and isoamyl alcohol
- Potassium thiocyanate (KCNS) solution: Dissolve 40 g of KCNS in water, and dilute solution to a volume of 100 ml
- Stannous chloride ($\text{SnCl}_2 \cdot 2\text{H}_2\text{O}$): Put 40 g of $\text{SnCl}_2 \cdot 2\text{H}_2\text{O}$ in 20 ml 6.5 N HCl. Add water to dissolve the salt and dilute the solution to volume of 100 ml.
- Wash solution: Take 10 ml of 6.5N HCl containg Fe^{3+} , 35ml of water, 1 ml of $\text{SnCl}_2 \cdot 2\text{H}_2\text{O}$, 2 ml of organic extractant, and shake for 2 minutes. Allow the phases to separate for 10-15 minutes. Use the aqueous phase of this solution as wash solution for removing interfering colour due to titanium and vanadium.
- Standard Mo solution, 100 ppm: Dissolve 0.150 g of molybdenum trioxide (MoO_3) in 100 ml of 0.1 N NaOH and dilute to little less than 1 litre. Make this solution slightly acidic with HCl and make the volume to 1 litre with distilled water. This solution should use for preparing working standards.

Procedure

Extraction

- Weigh 25 g soil in Erlenmeyer flask or Polyethylene bottle, and add 250 ml of acid ammonium oxalate solution. Shake for 16 hours. Then filter suspension through Whatman filter paper (Filter paper should initially washed with 6.5 N HCl) and discard the first 10-15 ml of filtrate.
- Transfer an aliquot of the filtrate 200ml to a porcelain evaporating dish.
- Evaporate the solution to dryness on a water bath.
- Ignite the residue with high temperature 500 °C in a muffle furnace for 5 hours to destroy oxalate and organic matter.

Estimation

- Cool the dish, and add 10ml of 6.5N HCl, and transfer the solution quantitatively to a 50 ml separating funnel for molybdenum determination.
- Dilute the HCl-molybdate to 45ml with distilled water, and add 1 g of citric acid and 2 to 3 ml of organic extractant, shake the mixture and allow 10-15 minutes for the phases to separate; discard the organic extractant which is the lower phase of the funnel.
- Again mix 1 ml of the KCNS solution with the sample solution, followed by 1 ml of the SnCl_2 solution.
- Add 1.5 ml of organic extractant to the coloured complex and shake for 2 minutes. Allow 10-15 minutes for phase separation, and then discard the aqueous phase, using fine tubing dropper.

[Handwritten signature]

[Handwritten signature]

[Handwritten signature]

[Handwritten signature]

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- Add 20 ml of wash solution and shake mixture for 2 minutes. Again allow 10-15 minutes for phase separation, and then discard the aqueous phase, using fine tubing dropper. Repeat this step, and, upon completing this shaking time, invert the funnel to facilitate washing out its tip with water.
- Dry the tip with the support of vacuum. Open the stopcock of the inverted funnel, and dry the stopcock bore with vacuum.
- Close the stopcock, shake the funnel quickly, and allow the phases to separate for about 15 minute by making funnel in upright position.
- The clear coloured supernatant solution is used for measurement (Mo-thiocyanate complex) at 470 nm.
- The work on blank and prepared standard solution (0.1, 0.4, 0.6, 0.8, 1 ppm) should be done as done above in the sample.

Calculation

$$M_o = C \times (250/25) \times (V_1/V)$$

C = concentration determined from standard curve

V_1 = Volume of the aliquot made upon following various treatments

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